SYALTEN® 30 mg

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Hyaluronic acid sodium salt for intra-articular use

PRODUCT DESCRIPTION:

SYALTEN 30 mg is a sterile, biodegradable, isotonic, injectable gel, for intra-articular use. **SYALTEN 30 mg** consists of medium chain (1.0-1.5 x 10⁶ Dalton) hyaluronic acid, obtained from *Streptococcus equi* bacteria, formulated to a concentration of 15 mg/ml in a physiologic buffer. **SYALTEN 30 mg** is characterised by viscoelastic properties, therefore allows to facilitate the normalisation of the viscosity of the synovial fluid present in the intra-articular cavity. Each box contains one syringe of **SYALTEN 30 mg** and a product leaflet. A set of two labels showing the batch number is contained in the box. One of these labels should be attached to the patient's file and the other should be given to the patient to ensure traceability.

COMPOSITION:

Sodium hyaluronate (15 mg/ml), sodium chloride, sodium dihydrogen phosphate dihydrate, dibasic sodium phosphate dodecahydrate, WFI grade water.

INDICATIONS:

SYALTEN 30 mg is a synovial fluid substitute which, thanks to its viscoelastic and lubricant properties, promotes the restoration of rheological conditions of the joints, altered in degenerative or post-traumatic conditions. The product, improving the characteristics of the synovial fluid, exerts a protective action of the joints and helps the improvement of joint function and the reduction of pain symptoms. **SYALTEN 30 mg** acts only at the joint where it is injected, without exerting any systemic action.

WARNINGS - PRECAUTIONS FOR USE:

SYALTEN 30 mg is suitable only for intra-articular injections and **must only be dispensed by a doctor who has received specific training on the intra-articular injection technique**. Before use, check the integrity of the syringe and the expiration date. Do not use needles other than those listed. The product should not be injected in the presence of a infected or severely inflamed joint. The infiltration must be avoided in the case of infections in place or inflammatory conditions of the skin in proximity of the injection. As no clinical experience is available for the use of Hyaluronan in children, treatment with **SYALTEN 30 mg** is not recommended in these cases. After the intra-articular injection it is advisable to recommend to the patient to avoid physical activities demanding stress for the articulation and resume normal activities after a few days.

SYALTEN 30 mg is a disposable product, the quality and sterility are guaranteed only if the syringe is sealed. Any residue must be discarded and not reused even after new sterilisation. Do not use the product if the package is already opened or damaged. The assembled syringe must be discarded immediately after use, regardless of whether or not the solution has been completely administered. After use, dispose according to applicable national practice.

INCOMPATIBILITIES:

There are incompatibilities between sodium hyaluronate and quaternary ammonium compounds, such as solutions of benzalkonium chloride. Contact between **SYALTEN 30 mg** and these substances should be therefore avoided.

SIDE EFFECTS:


There may be some temporary side reactions following injection of **SYALTEN 30 mg**, such as pain, stiffness, warmth, redness or swelling. These secondary manifestations may be relieved by applying ice on the treated articulation. Usually these effects disappear after a short time. If symptoms persist, consult a physician. Any other unwanted side effects associated with the injection of **SYALTEN 30 mg** should be reported to a doctor.

METHODS OF USE:


Remove any joint effusion before injecting **SYALTEN 30 mg**; for the removal of the effusion and the injection of **SYALTEN 30 mg** the same needle must be used. Remove the protective cap of the syringe, with particular attention to avoid contact with the opening. Firmly screw the needle, of diameter between 18 and 22 g, at the collar of Luer lock, following the instruction given below. Before injection, the site should be treated with appropriate disinfectant. Inject **SYALTEN 30 mg** adopting aseptic technique. Inject only into the joint cavity. It is recommended to perform an initial cycle of 3 weekly treatments, possibly followed by maintenance sessions, according to the medical prescription.

INSTRUCTIONS FOR ASSEMBLY OF THE SYRINGE NEEDLE:

- A. Carefully unscrew the cap of the tip of the syringe, being particularly careful to avoid contact with the opening.

- B. Gently grip the needle guard and mount the needle on the luer-lock mount, screwing it tight until a slight counter-pressure is felt in order to ensure an airtight grip and prevent leakage of the liquid during administration.

STORAGE:

Store **SYALTEN 30 mg** at 2–25°C (36–77°F) in a dry place in the original box. Protect from light, heat and frost. Keep out of reach of children.

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